USA/GB Technical data

5 Sizes:
XS (up to 36), S (37–38), M (39–41), L (42–44), XL (45–47)

- > available in knee-high and mid-calf versions
- > can be worn on the right or left leg

Material:
Outer sole: TPR (Thermoplastic Resin)
Inner sole: EVA
Inner material: Polyester
Shaft: Nylon

Cleaning instructions

To clean the shell and inner lining, please use lukewarm water (30° C possible) and mild soap if necessary. Wash the inside and dry immediately. Leave it in the open air to dry. Do not use a dryer with hot air.

Caution

Your doctor has prescribed you a stabilising orthosis. The design of the rocker sole warrants a physiologically correct gait by means of the »rolling-off« movement. Walking with the AirTraveler Walker is more difficult than walking with conventional footwear and its correct application is essential for the walker's functionality – please read the instructions.

- > Do not attempt to drive while wearing this walker.
- > Use extreme caution when climbing stairs, stepping up onto curbs and walking on uneven surfaces.
- > Walk more slowly than usual. You are not wearing a normal street shoe and more clearance is required for the wedge sole.

The application of this aid should be supervised by a physician or health-care specialist. In case of pain, swelling or skin irritation please call a doctor immediately!

Avoid inflating the air chambers too forcefully. Excessive pressure on the lower leg should be excluded!

Caution is required in patients with neuropathy to ensure that the walker does not restrict blood circulation. Patients should call their physician immediately if any tingling, burning, numbness or foul odour becomes noticeable.

If the product is improperly handled, its function and optimal safety can no longer be guaranteed.

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USA/GB Instructions for use

AirTraveler Walker

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1. Open the Velcro straps, remove the inner shoe of the walker and take out the hexagonal foam pad. (Figure 1)
2. Place the inner shoe around the leg and fasten it with the Velcro straps. (Figure 2)
3. To enter the walker safely, spread the medial and lateral splints using both hands. (Figure 3)
4. After entry, remove the white protective strips on the inside of the splint and fix the splint to the inner shoe. Now fasten the Velcro straps, starting with the straps close to the ankle. Ensure that the leg is comfortable in the walker. (Figure 4)
5. Remove the black pump on the back of the walker.
6. Now pump air through the valve on the back of the walker into the air chambers of the inner shoe. Use the transparent end of the pump (Figure 5). The air is distributed evenly. The degree of compression should ensure that neither the foot nor the ankle can move freely. The inner shoe of the walker should provide adequate protection for the foot and ankle.

Should impaired sensation, swelling or pain occur, the pressure should be reduced by releasing air from the chambers. Use the black end of the pump (Figure A). Should persistent pain occur, the treatment should be interrupted and a physician should be consulted.

Optional: Insertion of the foam pad (supplied as an accessory):

The four hexagonal pads (Figure B) are used to immobilise and pad the ankle, thus allowing individual fitting. Place pads in areas where friction needs to be reduced or excessive space needs to be filled.

Indications

- > Treatment of oedema
- > Acute management of ankle sprains
- > Injuries of the superficial tissues in the lower leg
- > Stress fractures
- > Fractures in the foot and ankle joint region
- > Fractures in the metatarsal region
- > Achilles tendon injuries
- > Injuries in the forefoot, midfoot and hindfoot region

Contraindication:

- > The AirTraveler Walker is not recommended when sensation in the foot and ankle region is impaired.

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