

MA Spacer (MIS Akin Spacer)

Toe separator

Ensure the outcome after minimally invasive Akin surgery without osteosynthesis material by using the new toe separator made from medical silicone.

The MA Spacer (MIS Akin Spacer) made from soft, viscoelastic silicone ensures the surgical outcome by fixing the corrected toes while adapting individually to the contours of the toes.

The biocompatible material is washable and can be disinfected, if required.

Indications

- > postoperative care of minimally invasive hallux valgus surgery (Akin osteotomy without osteosynthesis material)
- > conservative treatment of hallux valgus

New from DARCO



DARCO



Suggestion for postoperative care

one week after surgery
taping and use of the toe separator

for approx. 4–6 weeks
change the tape weekly and, while doing so,
clean and reuse the MA Spacer

from week 7
independent use of the toe separator by
the patient for approx. 4 weeks



The space between the first and second toes is preserved when narrow everyday shoes are worn again.



Purpose

- > ensures the surgical outcome after MIS Akin surgery with approx. 10° overcorrection
- > can also be used conservatively
- > corrects and separates toes that are close together
- > fixes and aligns toes that have undergone surgery / corrected toes in the desired position using tapes
- > provides protection from pressure pain

Features

- > biocompatible and non-slip
- > soft, smooth medical silicone for high wearing comfort
- > the material adapts individually to the contours of the toes
- > after wearing, it returns immediately to its original shape
- > easy to use
- > durable material
- > overcorrection of approx. 10°
- > washable and can be disinfected, if required, therefore more hygienic than other materials
- > cost-effective, as only 1 MA Spacer (MIS Akin Spacer) needed per patient
- > can be used on the left or right foot

Use

The MA Spacer (MIS Akin Spacer) is placed into the space between the first and second toes with the long part in a dorsal direction and is fixed with a tape strip.

SIZE

ART. NO.


Color: White | Universal size
Packaging unit: 20 pieces

MAS-20

SUPPLEMENTARY PRODUCT:

MECRON ReadyDressment Strips Elastic Stabilising Strips
(optionally available)

Ready-to-use, elastic pre-cut tapes are ideally suited for aligning toes that have undergone surgery / corrected toes and for fixing the MA Spacer.



Made in Germany

DARCO (Europe) GmbH | Gewerbegebiet 18 | 82399 Raisting | Germany
Phone +49 (0) 88 07.92 28-0 | Fax -22 | info@darco-europe.com | www.darco-europe.com

CE

190502