Body Armor[®] Pro Term

EN Foot Stump Orthosis
DE Fußstumpforthese
ES Ortesis para muñón de pie
FR Orthèse pour moignon de pied
IT Ortesi per moncone di piede
NL Voetstomporthese
DA Fodstumporthese
SV Fotstumpsortor
PL Orteza kikuta stopy
SL Nogavična ortoza
HR Ortoza za bataljak stopala
PT Órtese de coto do pé

DARCO

1.

2.

3.

4.

5.

6.

Body Armor® Pro Term | Foot Stump Orthosis

Important notes

If there is a serious deterioration in your state of health as a result of using the medical device, please report this to your specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, the dealer will forward this immediately. The reports are agreed with and processed by the competent authorities.

Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This aid assists recovery and should be used exclusively under the instruction of and in consultation with your doctor or a healthcare professional. DARCO (Europe) GmbH does not assume any responsibility. Do not wear the product with injured, irritated, damaged or intact skin. You should use bandages or socks, for example.

We recommend that you do not drive a vehicle while you are wearing this aid. Use extreme caution when climbing stairs (e.g. also stepping up onto curbs) and walking on uneven surfaces. Walk more slowly than usual.

This aid is intended for single-patient use only.

Purpose

Foot Stump Orthosis as an interim treatment for stabilisation and offloading with simultaneous mobilisation of the patient.

Indications / Contraindications

Indications

- > Diabetic foot syndrome
- > PAOD (peripheral arterial occlusive disease)
- > Trauma
- > Postoperatively for treatment following Chopart amputation, Lisfranc amputation

Contraindications

No contraindications are known.

Design features

Stable, ankle-spanning shaft with wedge-shaped insole for offloading of the stump. The integrated and adjustable air cushions enable individual adjustment. The product includes a ventral support strap with adjustable closure system for good stability.

Technical data

2 sizes

M = Stump length 16 cm

L = Stump length 19 cm

Universal fit for right and left.

Colour

Black-Blue

Material

Outsole walking surface: TPR

Outsole body: EVA

Plastic shell: Polyurethane

Fabric liner: Polyester

Buckles: Polypropylene

Insole: EVA, PE

- > **Photo 1:** Open the two fasteners on the orthosis and remove the blue inner shoe. Open the hook and loop fasteners of the inner shoe and place this around the leg. Fasten the inlay using the hook and loop fasteners. First close the hook and loop strap over the instep, then the lower hook and loop strap as numbered in the photo.
- > **Photo 2:** To insert the foot securely into the walker, spread the medial and lateral parts of the shell with both hands. Now place both insoles into the walker and step onto the two insoles with the blue inner shoe. The insoles can be adjusted individually by qualified staff.
- > **Photo 3:** After inserting the foot, pass the straps into the buckles and adjust them to the desired pressure. Fasten the buckles.
- > **Photo 4:** Now pump air through the round pump on the front of the walker into the air chambers of the blue inner shoe. The air is distributed evenly. The compression achieves additional stability of the ankle joint.
- > **Photo 5:** To reduce the volume of the air chambers, release the air using the silver valve.

Note

If impaired sensation, swelling or pain occurs, the pressure should be reduced by releasing air from the chambers. Should persistent pain occur, the treatment should be interrupted and a doctor should be consulted.

Combination with other products

No specifications

Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling/cleaning. The medical device is designed to be used for the usual duration of use prescribed by a doctor without showing a reduction in safety or performance. After the period of use prescribed by your doctor, check the function of the product together with your doctor or specialist dealer.

Maintenance

The product does not require maintenance.

Cleaning instructions

Hand wash

If the shoe is slightly soiled, you can wipe it with a damp cloth or clean it with a soft brush. To clean the shell and liner, use lukewarm water (30° is possible) and mild soap or mild detergent if necessary. Rinse and dry the liner promptly. Do not dry this in a dryer or on a radiator. Let it dry in the air at room temperature.

Disposal

Dispose of after use. Please dispose of the packaging in an environmentally friendly manner. Dispose of this product after use in accordance with applicable local laws and regulations.

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