**PegAssist® Insole**

Off-loading insoles

The PegAssist® Insole is a system that effectively off-loads the plantar aspect of the foot. The easy removal of single pegs allows for an individual customization in post-operative care or acute wound treatment. The thin cardboard cover stabilizes the off-loaded area and prevents the surrounding pegs from shifting. The pressure is reduced by up to 60%.

**Recommendations for use**

1. Remove the black original insole from the shoe (fig. 1).
2. Insert in return the PegAssist® Insole into the shoe with the blue side facing upwards (and the pegs facing downwards) (fig. 2).
3. Place a sterile film over the wound or pressure ulcer (fig. 3) and draw a circle on the film around the affected area with the included marker pen (fig. 4).
4. The patient should now stand up on the insole in the shoe so that the marked area is transferred onto the blue upper layer of the insole (fig. 5).
5. Carefully remove the pegs beneath the marked area where the off-loading will take place (fig. 6).
6. Apply the complete surface of the cardboard cover (included in the pack) to the insole with the pegs. The cardboard has an adhesive backing which stabilizes the pegs surrounding the off-loaded area and ensures non shifting (fig. 7).
7. Reinsert the adjusted PegAssist® insole (with the cardboard cover facing downwards) into the shoe (fig. 8).
8. The system is ready (fig. 9).

**Caution**

It is recommended to wear bandages or socks. If the product is improperly handled, its function and optimal safety can no longer be guaranteed. Signs of fatigue can occur in the cover layer when there is a large gap. This aid is intended for single-patient use only.

**Indication**

Selective off-loading for:

- > plantar lesions
- > plantar ulcerations
- > pressure points

**Contraindication**

- > open fractures

**Technical data****Size:**

- > available in seven sizes
- > right and left fit

The insole is available with two different colors of pegs (blue or black – depending on the model) in all required sizes for the following foot and ankle products:

- > OrthoWedge – Off-loading Shoe
- > OrthoWedge Light – Off-loading Shoe
- > Relief Dual® – Off-loading Shoe
- > HeelWedge® – Off-loading Shoe
- > MedSurg, MedSurg Pro® – Healing Shoe
- > SlimLine® – Cast Boot
- > Relief Insert® Walker – Lower leg / foot orthosis

**Material:**

The 18 mm thick insole consists of 3 layers: a combination of EVA (blue or black), Plastazote® and a top layer of Poron®.

**Maintenance:**

The product requires no maintenance.

**Cleaning instructions:**

The insoles can be cleaned with warm water and soap, or with commercial disinfectants.

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200213

Made in China

Size:	REF:	Size:	REF:	Size:	REF:
XL	PAS4	XL	PAS4	XL	PAS4
L	PAS3	L	PAS3	L	PAS3
M	PAS2	M	PAS2	M	PAS2
S	PAS1	S	PAS1	S	PAS1
XS	PAS0	XS	PAS0	XS	PAS0

for / für OrthoWedge, OrthoWedge Light, HeelWedge®, SlimLine®  
for / für MedSurg, MedSurg Pro®, Relief Dual®, Relief Insert® Walker

The DARCO PegAssist® System features a multi-indication removable-peg chassis that effectively off-loads the plantar aspect of the foot after surgery or when wounds are present.

> Removable pegs allow for localized off-loading of wounds and ulcerations of the foot.

> Pressure is reduced by as much as 60% while allowing patient to remain ambulatory.

> Included stabilizer board prevents adjacent pegs from collapsing.

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