


Bungee Wrist Splint – Comfortable and Easy to Use

New at DARCO



- > mild sprains and strains
- > tendonitis and carpal tunnel syndrome
- > injuries and tenosynovitis
- > interim provision after cast removal
- > carpoptosis, wrist drop
- > rheumatoid arthritis, stabilizing effect



DARCO



Bungee Wrist Splint

Eight elastic bungees attached to a single Velcro® strap provide:

- > compression
- > contoured fit and
- > easy application

Achieve complete immobilization the wrist with the aluminium palmar stay. The aggressive palmar styling allows for full function of all fingers. The brace is fully washable and comes with a laundry bag.

Fabric

Comfort equals compliance. Tricot padded lining, perforated sueded exterior and nylon trim all make the Bungee Wrist Splint comfortable to wear. A dorsal stockinette gusset prevents lace contact with the skin and aids in one-handed application.

Fastening

Easy to use bungee strap allows the wearer to apply, remove and make quick adjustments without assistance.

Item-No.	Colour	Size	Wrist Circ.	Length
left				
BW0-L	Black	XS	13.34 – 15.88 cm	20.3 cm
BW1-L	Black	S	15.88 – 18.41 cm	20.3 cm
BW2-L	Black	M	18.41 – 20.32 cm	20.3 cm
BW3-L	Black	L	20.32 – 22.86 cm	20.3 cm
BW4-L	Black	XL	22.86 – 25.40 cm	20.3 cm
right				
BW0-R	Black	XS	13.34 – 15.88 cm	20.3 cm
BW1-R	Black	S	15.88 – 18.41 cm	20.3 cm
BW2-R	Black	M	18.41 – 20.32 cm	20.3 cm
BW3-R	Black	L	20.32 – 22.86 cm	20.3 cm
BW4-R	Black	XL	22.86 – 25.40 cm	20.3 cm