

DARCO TAS[®] Toe Alignment Splint for postoperative treatment – also ideal after minimal invasive interventions



Modular product elements



High wearing comfort, even under socks



The DARCO TAS[®] Toe Alignment Splint is a cost-effective solution that provides superior post operative toe alignment after Hallux Valgus ⁽¹⁾, Hammer Toe ⁽²⁾ and Tailor's Bunion ⁽³⁾ procedures.



DARCO

TAS[®] customizes to each patients needs!



The TAS[®] is an effective and inexpensive way of maintaining correct alignment of the toes during the healing process.

A thin elastic metatarsal band reduces slippage and soft toe straps hold the toes in the correct position to maintain alignment following surgery to correct hallux valgus, hammer toes or Taylor's bunions. May also be used in the treatment of capsulitis plantar to the m.p. joints caused by mild subluxations of the joint. The TAS[®] is a well-established standard treatment following forefoot surgeries.

Features and Benefits

- > Low profile metatarsal band fits easily into most footwear and reduces slippage.
- > Soft toe straps will not irritate skin.
- > Special T-Strap gently guides toe(s) into proper alignment for post-procedural healing.
- > One sizes fits all.
- > Latex Free

Item-No.	Colour	Size (UK)
----------	--------	-----------

TAS [®] (Toe Alignment Splint)	White	Universal
---	-------	-----------

right and left fit

