

PegAssist® Insole

For selective off-loading and acute treatment



- > removable pegs for effective off-loading
- > shock-absorbing effect
- > available for the shoe models shown on reverse side
- > right and left fit
- > variety for OrthoWedge Light



DARCO

Pegs can be removed from specific areas of the insole in order to redistribute pressure. The cardboard cover then stabilizes the off-loading area.






The PegAssist® Insole is available for all products listed here.

The 18 mm thick insole is layered with Poron®, Plastizote® and EVA to provide the perfect combination of shock absorption and moldability.

After removing the pegs from specific areas of the insole, the adhesive-backed cardboard cover is then applied to the complete surface of the insole. The cardboard cover acts to stabilize the remaining pegs.



	SIZE (D/EU)	ITEM-NO.					
	for MedSurg /-Pro®, Relief Dual®			for OrthoWedge, HeelWedge®, A.P.B.®, SlimLine®			
	Women						
	WS	34.0 – 36.5		PAS-W1	XS	33.5 – 37.0 PAS0	
	WM	37.0 – 38.5		PAS-W2	S	37.5 – 39.5 PAS1	
	WL	39.0 – 41.0		PAS-W3	M	40.0 – 41.5 PAS2	
					L	42.0 – 44.0 PAS3	
					XL	44.5 – 47.0 PAS4	
	Men						
	MS	39.0 – 41.0		PAS-M1		for OrthoWedge Light	
	MM	41.5 – 43.0		PAS-M2		XS	33.5 – 37.0 PAS-STL0
ML	43.5 – 45.0	PAS-M3	S	37.5 – 39.5 PAS-STL1			
MXL	45.5 – 47.0	PAS-M4	M	40.0 – 41.5 PAS-STL2			
			L	42.0 – 44.0 PAS-STL3			
			XL	44.5 – 47.0 PAS-STL4			

