

Keep the balance: TwinShoe balance shoe by DARCO.



After foot surgery, most surgeons employ a post-operative shoe to enable the patient to become mobile as quickly as possible. The consequence of wearing such a shoe is a height discrepancy between the two legs which results in discomfort, imbalance and possible danger when walking.

DARCO's TwinShoe balance shoe lessens the height discrepancy and makes walking much easier and safer.



DARCO



TwinShoe balance shoe

The TwinShoe is a matching and balancing shoe designed to equalize the different leg lengths that occur when wearing most post operative shoes. The TwinShoe is a perfect complement to our well-known products such as the OrthoWedge and OrthoWedge Light forefoot offloading shoes, HeelWedge® heel offloading shoe, Body Armor® Walker II cast replacement and the WCS™ long-term healing shoe.

The TwinShoe makes it no longer necessary to spend time applying and removing complicated height adjustments. Simply slip it on like any other shoe and go!

The TwinShoe is designed for a universal left / right fit and is the matching shoe for the following DARCO products:

Item-No.	Size (D/EU)	Colour
TS0B	XS 33.5 – 37.0	Black
TS1B	S 37.5 – 39.5	Black
TS2B	M 40.0 – 41.5	Black
TS3B	L 42.0 – 44.0	Black
TS4B	XL 44.5 – 47.0	Black

